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



CADUCIDAD: 12 meses. Conservación a: 2°C-8°C

Ed. 020217

INgezim PRRS DR is based on a double recognition ELISA technique, which uses the PRRSV's N protein both, as antigen and as conjugate.

TECHNICAL BASIS OF THE KIT

- Plates are coated with recombinant European PRRSV's N protein. Conjugate and serum samples are added to the wells and incubated.
- The recombinant European N protein HRPO conjugated will bind to the antibodies present in the sample if the serum is positive. This binding is detected by adding a specific substrate.

APPLICATION

Due to its characteristics of sensitivity and precocity, the assay can be a useful tool to control the serological status of new animals and to monitor negative herds. Moreover, the assay is able to detect late humoral responses. The kit has been developed using the European PRRSV N protein as antigen. Nevertheless, due to the homology existing between both kinds of strains, the kit must not be used to differentiate between American and European strains.

INTERPRETATION OF THE RESULTS


There is one Cut Off value for the results interpretation. Samples with an OD higher than the Cut Off Value are positive and samples with an OD lower than the Cut Off value are negative.

ANALYTICAL SENSITIVITY

STUDY 1: 6 animals experimentally infected and bled at different days post challenge until 60 days post infection; and 69 pigs experimentally infected with European strain of PRRSV and bled at days 3, 7, 14, 21 and 28 post infection.

STUDY 2: 38 four weeks old piglets experimentally infected with 5 different strains of genotype I and bled till day 28 p.i. (CRESA, Centre de Recerca en Sanitat Animal, UAB-IRTA, Bellaterra, Spain).


STUDY 3: Sera of animals experimentally infected with a European strain and bled at different days post infection.

The results obtained indicated that the assay is able to detect specific antibodies to protein N at day 7 post infection in 88.4% of cases in study 1, 76% in study 2 and 93.4% in study 3.

DETECTION OF LATE HUMORAL RESPONSE.

16 animals experimentally infected with 2 different strains of genotype 1 and bled till day 84 post infection. The assay maintains the ability of detection at day 84 p.i.


J of Virol Methods 181 (2012) 109-113 ; JVDI 2012 Mar; 24(2):344-8

COMPARATIVE WITH OTHER CONVENTIONAL KITS

STUDY 1: 48 positive sera by PCR. The assay detected 13% more positive than HerdCheck 2xR.

STUDY 2: 406 sera from 8 positive herds. The assay was able to detect from 14% to 50% more positive samples than the rest of ELISAs.

STUDY 3: Sera of animals infected with 8 different European strains (serotypes I and III). The assay detected 10% more positive than HerdCheck 3x and 4% more than HerdCheck 2xR. (Piwet, Pulawy).

STUDY 4: Sera of animals infected with 6 different European strains (Lelystad, 07V063, 08V204, LENA, 08V194, OLOT). The assay detected specific antibodies at day 7 p.i. except for strain 08V204 (day 10 p.i. and OLOT (day 14 p.i.).

SPECIFICITY

1002 sera from negative herds without clinical symptoms (819 Spanish and 183 Canadian). The specificity was **99%**.

RING TEST EVALUATION

In 2010 and 2011, the assay was evaluated in the "International PTS PRRSV antibody detection in serum" organized by GD B.V. In both years, the assay resulted to be more sensitive than the rest of assays available in the market for European strains.

REFERENCES

- J. of Virol. Methods 181 (2012) 109-113
- JVDI 2012 Mar; 24(2):344-8

IT-73840
IT-73780
ISO 14001:2015
9191.INGE
ISO 9001:2015
9175.ING2

Shelf life: 12 months. Store at: 2°C-8°C

Ed.020217