

## INgezim® TGEV 2.0

R.11.TGE.K3

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		INgezim DIFFERENTIAL CORONAVIRUS				
		PRVC	TGEV	NEGATIVE	DOUBTFUL	TOTAL
COMPETITI ON ELISA	PRCV	498			5	503
	TGEV	1	31		1	33
	NEGATIVE	10		112		123
	DOUBTFUL	6				6
	TOTAL	516	31	112	6	665

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

## INgezim® TGEV 2.0

R.11.TGE.K3

**INgezim® TGEV** is based on an enzymatic immunoassay (Capture Blocking ELISA) which uses as conjugate a monoclonal antibody (MAb) specific to site Ac of TGEV S protein (specific epitope of TGEV)<sup>1</sup>. An additional MAb is used as a capture reagent of the S protein.

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with recombinant S protein of TGEV which is captured by a specific monoclonal antibody. Each serum sample must be added to the well.
2. If the serum sample contains specific antibodies to the Transmissible Gastroenteritis Virus (TGEV), they will bind to all of the epitopes present.
3. After an incubation period, a specific TGEV Mab-HRPO (Conjugate B) is added to the wells. They will only bind to free epitopes. The presence or absence of the labeled Mab can be detected by adding a substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

### APPLICATION

Detection of Transmissible Gastroenteritis Virus (TGEV) specific antibodies in porcine serum samples. There is no interference with any other porcine coronaviruses.

### INTERPRETATION OF THE RESULTS

Two cut off are used for the results interpretation: positive & negative to TGEV.

### VALIDATION

#### EXTERNAL EVALUATION

##### Comparison with seroneutralization (SN).

- A study with 411 samples, divided into 4 groups, was made:
- 94 samples from 16 herds serologically negative to PRCV and TGEV.
  - 124 samples from 21 herds serologically positive to PRCV.
  - 122 samples from 17 herds serologically positive to TGEV
  - 71 samples from 22 diagnostic cases

The results obtained indicated **94% sensitivity** and **98.2% specificity**.

##### Analytical sensitivity<sup>2</sup>.

Samples of 30 piglets taken off at days 0, 4, 7, 11, 14 & 18 p.i. were analyzed. Results obtained indicate that the assay is able to detect antibodies specific of TGEV since day 14p.i. depending on the animal studied.

#### INTERNAL EVALUATION

##### Comparison with an ELISA that uses complete TGEV virus as antigen instead of recombinant S protein.

A set of 665 sera were analyzed and the results obtained in both assays were compared.

		INgezim DIFFERENTIAL CORONAVIRUS				
		PRVC	TGEV	NEGATIVE	DOUBTFUL	TOTAL
COMPETITIVE ON ELISA	PRCV	498			5	503
	TGEV	1	31		1	33
	NEGATIVE	10		112		123
	DOUBTFUL	6				6
	TOTAL	516	31	112	6	665

Those samples with discrepant results showed OD values close to the cut off when in the ELISA that used complete virus as antigen. The use of recombinant protein as antigen allows the detection of specific antibodies to TGEV, avoiding false positive animals.

<sup>1</sup>Developed in the Laboratory of Dr. Luis Enjuanes (CNB, Madrid). Carlos M. Sánchez et al. Antigenic Homology among Coronavirus related Transmissible Gastroenteritis Virus. Virology 174, 410-417.

<sup>2</sup>Chai et al. BMC Veterinary Research 2014, 10:75

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells.
- Vials with Positive Control for TGEV
- Vials with Negative Control for Porcine Coronavirus
- Vials with Peroxidase TGEV specific Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA  
REGISTRATION NUMBER 10477 RD

**SHELF LIFE: 18 months**  
Stored at 2-8°C

Ed.010818