INgene q Brucella Sondas Typing

R.10.BRM.K.5TC/Q

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PCR Real-time. 100 Reacciones / 100 Reactions

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Table of Contents (ENGLISH)

1.	INTRODUCTION.....	8
2.	KIT COMPOSITION.....	8
3.	GUIDELINES FOR THE CORRECT CONSERVATION OF THE KIT COMPONENTS.....	8
4.	MATERIALS REQUIRED BUT NOT SUPPLIED.....	8
5.	PRECAUTIONS TO AVOID CONTAMINATION.....	8
6.	SAMPLES COLLECTION AND TRANSPORT.....	9
7.	DNA EXTRACTION.....	9
8.	AMPLIFICATION OF THE GENETIC MATERIAL.....	10
9.	ANALYSIS AND INTERPRETATION OF RESULTS.....	11
10.	TROUBLESHOOTING.....	12
11.	QUALITY CONTROL.....	13
12.	TECHNICAL ASSISTANCE.....	13

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1. INTRODUCTION

INgene q Brucella Typing is a real-time PCR kit suited for the detection and typing of DNA of *Brucella* (BRM) in biological samples in a simple way, with high levels of sensitivity and specificity.

The kit contains all necessary reagents and enzymes in two mixes. Once mixed, they are at the required concentrations to detect and type (and optionally quantify) BRM, simply by adding the DNA of the problem sample. The kit amplifies several parts of the genome, which differ according to species. This kit makes use of three TaqMan-MGB-NFQ probes marked with FAM, VIC and TEXAS RED.

The kit includes a BRM positive control for detection. The kit does not include a standard positive control with a known number of BRM copies (ask our commercial department for standard control). This control allows generating a standard curve which links the number of pathogen copies with the cycle threshold (Ct) value.

The kit has been tested for ROCHE LC480. It is possible to work with other thermocyclers from other commercial brands, but Ct values may vary among them.

CAUTION: follow specific software of each system.

2. KIT COMPOSITION

COMPONENTS	Nº VIALS	VOLUMEN/VIAL
Mixture A (BRM specific primers and probes)	2	600 µl
Mixture B (enzyme mixture)	2	600µl
Positive Control A1 – BRM Amplification control	1	60 µl

3. GUIDELINES FOR THE CORRECT CONSERVATION OF THE KIT COMPONENTS

On receiving the kit, keep it at -20°C until use, except mixture B, which must be kept at 4°C. The components of the kit are stable for 1 year from the manufacturing date (see expiry date on packaging).

4. MATERIALS REQUIRED BUT NOT SUPPLIED

- Dust-free disposable gloves (without talcum powder).
- Microcentrifuge.
- Tube shaker.
- Real-time thermocycler.
- Micropipettes (0.5-1000 µl).
- Sterile pipette tips (with filter).
- Sterile DNAses / RNAses free water.

5. PRECAUTIONS TO AVOID CONTAMINATION

The following points should be read with care:

- Disposable items must be DNase and RNase-free.
- Use DNase and RNase-free distilled, autoclaved water (25 min., 120°C).
- Use sterile filtered tips.
- Ensure good homogenisation of the kit components once defrosted.
- Maintain the A and B mixes stored in ice at all times.
- Repeated cycles of freezing and defrosting may reduce the sensitivity of the reagents. Protect them from exposure to light until use.

To avoid contaminations giving rise to false positives it is important to:

- Physically separate the positive PCR control from the remaining reagents of the kit.
- To carry out any handling of samples to be tested in a different location/room from the one where the amplified products are being analyzed.
- Add the positive PCR control in a different location/room from the one where the mix is added and where the samples to be tested are being handled.

6. SAMPLES COLLECTION AND TRANSPORT

The assay could be performed using any biological sample considered of interest by the clinician. Our technique has been standardised and contrasted with a considerable number of liquid or colony culture, blood, sera, lymph node and milk. Samples must be cooled from the moment they are obtained until they are processed. Processing needs to occur within 24 hours after obtaining the samples. In case a longer time is envisaged, we recommend freezing the samples.

7. DNA EXTRACTION

Any extraction procedure should be used which yields good-quality of the extracted material.

In case you have no optimized extraction procedures, INGENASA could provide extraction protocol and buffer under request.

Ingenasa has tested the following methods for the DNA extraction:

- Method of Chomczynski and Sacchi (Phenol extraction).
- Automatic magnetic methods (Mag MAX™ Total NA isolation kit, Ref AM1840).
- Colum purification methods: QIAamp® Viral RNA Mini Kit (QIAGEN), High Pure PCR Template Preparation kit (ROCHE).

There are several commercial houses that provide high quality extraction kits that could be used following the manufacturer instructions, like: Life-technologies, ROCHE, Qiagen, Marchery-Nalgene, others.

SUGGESTED PROCEDURES

a) DNA purified by " NucleoSpin Tissue" , Marchery-Nalgen. Follow commercial house instructions.

b) Bacteria that have been grown in solid media (agar plates) must be used as starting material. Resuspend a bacterial colony with a seeding inoculation loop in 200 µl of sterile ultra-pure water. Incubate the suspension in a bath at 100°C for 10 min. Cool in crushed ice for 5 min. Use 1 µl of this suspension in the amplification mixture for PCR.

c) Alternatively, starting from a culture in liquid media is possible (10ml TBS culture growing with aeration at 37°C during 24-48 hours). Centrifuge (12,000 rpm, 5 min) 1 ml of the culture in an Eppendorf-type tube, remove the supernatant and resuspend the pellet again by means of stirring with the bacteria in 1 ml of sterile ultra-pure water. Centrifuge again, remove the supernatant and resuspend the pellet again in 1 ml of sterile ultra-pure water. Incubate the suspension in a bath at 100°C for 10 min. Cool in crushed ice for 5 min. Use 1 µl of this suspension in the amplification mixture for PCR.

CAUTION: If amplification is not going to be done immediately, preserve the DNA at -20 ° C.

8. AMPLIFICATION OF THE GENETIC MATERIAL

REQUIRED MATERIALS

- Crushed ice.
- DNA extracted from the samples.
- Mixtures A and B – **KEEP IN CRUSHED ICE AT ALL TIMES.**
- Positive amplification control A1 (BRM).
- DNase-RNase-free water.

PROCEDURE

1. Prepare and identify as many tubes for the amplification as samples to be processed, adding an additional tube for the positive amplification control, and another one for the negative control.
2. Take mixtures A and B out, keeping them in crushed ice, until complete thawing of mixture A. Make sure that they are correctly homogenised before taking out the required volume for the assay.
3. Prepare an appropriate amount of amplification mixture for the number of samples to be processed. The volume of each reagent to be mixed for each of the samples is:

	Mixture A	Mixture B	Final Master Mix Volume
Per sample	10 µl	10 µl	20 µl
For 10 samples	100 µl	100 µl	200 µl

The tube used for mixing should be kept in crushed ice at all times. Likewise, it is recommendable to prepare an excess amount of mixture (calculate an extra 10% for all reagents) in order to compensate for possible losses of volume during pipetting.

4. Once mixture is prepared, homogenise correctly. Place the tubes previously labelled in crushed ice and add 20 µl of the mixture prepared in this way to each tube.

Add 1 µl of previously extracted DNA samples to each tube, 1 µl of positive control A1 (amplification control for BRU), to the corresponding tube and 1 µl of water to the tube labelled as negative control. Carefully mix the contents of each tube and ensure that all the liquid has been deposited at the bottom of the tube. If not, tubes may be centrifuged lightly until this occurs.

5. Set the thermocycler to the following conditions:

	Temperature (°C)	Time	Cycles
Denaturation	95	10 min	1
Amplification	95	15 s	40
	60 *	60 s	

6. Reading the fluorescence takes place during the elongation step (marked in the table)*, and the channels through which the fluorescence data are collected are detailed in table:

	Reporter	Quencher
BRM	FAM	NFQ
	VIC	NFQ
	TEXAS RED	NFQ
Passive reference	None	

9. ANALYSIS AND INTERPRETATION OF RESULTS

Analysis of results

Each problem sample analyzed yields data regarding the fluorescence originated from channels FAM, VIC and TEXAS RED. Automatic threshold analysis (threshold and baseline) is recommended.

If this option is not viable, the analysis can be carried out manually following the thermocycler instructions.

FAM, VIC and TEXAS RED channels: pathogen detection

A positive result during amplification implies a typical fluorescence curve with a Ct value <40.

Samples with Ct's comprised between 38 and 40 should be regarded as doubtful. In that case we recommend:

- Repeating the analysis.
- Sequencing the amplified product.

Qualitative analysis:

Possible results are summarised in the table:

Table: Results validation

Probe	B.suis	B. melitensis	B.abortus	B.ovis
FAM	+	-	-	-
VIC	+	+	-	+
TEXAS RED	+	+	+	-

Case 1: there is amplification signal in all the possible channels.

Result: **the sample contains *Brucella suis*.**

Case 2: there is signal in channels VIC and TEXAS RED, but not in FAM.

Result: **the sample contains *Brucella melitensis*.**

Case 3: there is signal in channel TEXAS RED only.

Result: **contains *Brucella abortus*.**

Case 4: there is signal in channel VIC only.

Result: **contains *Brucella ovis*.**

10. TROUBLESHOOTING

Detection of fluorescence signal on any channel in the negative extraction control.

Possible cause	Solution
Contamination during extraction process	<ul style="list-style-type: none"> Repeat the DNA extraction process and PCR using new reagents. Ensure that the workplace, equipment and all instruments used for performing the PCR are decontaminated, i.e. free from nucleic acids.

Detection of fluorescence signal on any channel in the negative PCR control.

Possible cause	Solution
Contamination during PCR preparation.	<ul style="list-style-type: none"> Repeat PCR reaction with new reagents and, if possible, various replicates of each sample. If tubes are being used, close each one immediately after adding the sample. As a strict rule, the positive control must be added last and in a separate physical location from where the problem samples are added. Ensure that the workplace, equipment and all instruments used for performing the PCR are decontaminated, i.e. free of nucleic acids.

Absence of fluorescence signal on any channel in the positive PCR controls.

Possible cause	Solution
Wrong fluorochrome channels selected	<ul style="list-style-type: none"> Check that you select the FAM, VIC and TEXAS RED channels for all samples and controls to be analyzed.

11. QUALITY CONTROL

Each production batch of this kit has been checked under our ISO-certified Quality Management System.

12. TECHNICAL ASSISTANCE

For further information concerning the assay and its performance, please contact:

India Contact:
Life Technologies (India) Pvt. Ltd.
Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com